



EUROfusion PWIE SP-B.4

Reference Coatings for iTER and DEMO

IST activities for 2021: RBS and NRA characterization of selected Be and W reference samples - plans and capabilities

**Rodrigo Mateus
Norberto Catarino
Rui Silva
Eduardo Alves**




ipfn

**INSTITUTO DE PLASMAS
E FUSÃO NUCLEAR**



This work has been carried out within the framework of the EUROfusion Consortium and has received funding from the European Union's Horizon 2020 research and innovation programme under grant agreement number 633053. The views and opinions expressed herein do not necessarily reflect those of the European Commission.

Laboratory layout




- JET chamber - EBS, RBS, NRA, PIXE, from small samples to tiles**
- RBS/C, ERD chamber - EBS, RBS, NRA, Be-based, W-based, small samples, good adhesion**


Experimental Setup - 1




JET line





Experimental chamber



Detector geometry



RBS/C, ERD chamber



Detectors

- 1 Si(Li) detector / for NRA
 - Active area / 200 mm²
 - Active thickness / 2 mm
 - Absorber foil
 - 71900 x 10¹⁵ at/cm² Mylar (C₁₀O₄H₈)
 - Solid angle / 35 msr
 - Reaction angle / +140°
- 3 Si pin-photodiode / for RBS, EBS
 - Active thickness / 0.3 mm
 - Solid angles / 8.8 msr, 11.7 msr, 11.7 msr
 - Scattering angles / -140°, -165°, +165°

Case studies: Deuterium retention in Be-based films (IAP)

Table 1. List of Be-based coatings deposited on W plates at CT and up to 400 °C.

Reference	Nominal composition	Nominal thickness (μm)	Substrate temperature (°C)
60190516_1_2	Be+D(10%)	5	Chamber Temp. (CT)
60190520_1_2	Be+D(5%)	5	Chamber Temp. (CT)
60190523_1_2	Be+D(5%)	5	400
60190527_1_2	Be+D(10%)	5	200
60190528_1_2	Be+D(5%)	5	100
60190607_1_2	Be+D(10%)	5	JET-like pulses

WP PFC 2019





Table 2. Selected Be-based coatings.

Reference	Nominal comp. (at.%)	Nominal thickness. (μm)	Substrate	Substrate temp. (°C)
60200211_3_3	Be+D(5)+Ne(2.5)	5	W	Not heated
60200129_1_1	Be+D(5)+Ne(2.5)	5	W	100
60200130_1_2	Be+D(5)+Ne(2.5)	5	W	400

WP PFC 2020




Case study: Helium load on W-O coatings (ENEA)




Coating	Description	Nominal thickness (μm)
WFW536	compact W-O	>1.0
WFW549	porous W-O	>1.0

Ion beam	Ion implanted fluence
150 keV $^4\text{He}^+$	2×10^{17} ion/cm ²
100 keV $^4\text{He}^+$	2×10^{17} ion/cm ²
50 keV $^4\text{He}^+$	1×10^{17} ion/cm ²




Ion implantation



SRIM simulation for W:O (37:62)



TOF-ERDA (RBI) compact W-O



EBS / RBS porous W-O



Case study: Deuterium load on W-O coatings (ENEA)




Coating	Description	Nominal thickness (μm)
WFW600	porous W-O	>1.0
WFW624	porous W-N-O	>1.0


Ion beam	Ion implanted fluence
150 keV D_2^+	2×10^{17} ion/cm ²
100 keV D_2^+	2×10^{17} ion/cm ²
50 keV D_2^+	1×10^{17} ion/cm ²

WFW600

NRA




SIMS (VTT)




compatible results:
EBS/RBS/NRA (IST)
TOF-ERDA (RBI)
SIMS (VTT)

WFW624

NRA



SIMS (VTT)



Within the PWIE of the European Fusion programme IST will contribute with Ion Beam Studies of Be-based and W-based reference samples.

Besides EBS, RBS and NRA, IBA analysis at IST may easily include PIXE and PIGE.

Additional ion implantation experiments (^1H , ^2H , ^4He ...) are commonly performed.

