



EUROfusion

Work Package **PWIE**

Subproject E: PWI with Be, T and neutrons focus on JET post-mortem analysis and its interpretation

Kickoff Meeting : IST activities for 2022

Norberto Catarino

Rui Silva

Eduardo Alves



ipfn

INSTITUTO DE PLASMAS
E FUSÃO NUCLEAR



This work has been carried out within the framework of the EUROfusion Consortium and has received funding from the European Union's Horizon 2020 research and innovation programme under grant agreement number 633053. The views and opinions expressed herein do not necessarily reflect those of the European Commission.


SP E.2 Comparison of hydrogenic retention quantification by different techniques and fuel removal assessment

Del 5: Characterization of JET divertor tiles 0 and 1 using ion beam analysis (RBS, NRA).

SP E.3 Post-mortem analysis of PFC and other objects in JET


Del 6: Characterization of JET plasma facing and diagnostics components using ion beam analysis (RBS, NRA)

1. Laboratory Layout



Experimental Setup


JET line



Experimental chamber



Detector geometry



- **D profiling NRA**

Nuclear Reaction Analysis - $\text{D}({}^3\text{He}, \text{p}){}^4\text{He}$

- **Be and C Profiling**

Nuclear Reaction Analysis - ${}^9\text{Be}({}^3\text{He}, \text{p}){}^{11}\text{B}$ ${}^{12}\text{C}({}^3\text{He}, \text{p}){}^{14}\text{N}$



**Elemental Profiling: Rutherford Backscattering Spectrometry (RBS)
and Elastic Backscattering Spectrometry (EBS)**

- **Trace impurities (metals)**

PIXE (X-ray emission)

- **Computational analysis: WiNDF**

Case studies: JET results




(foto during 2014 shutdown)




Case studies: JET results (Tile 0)



14N HFGC LH ILW-3



14N HFGC LH [F630 00 002]





(foto during 2014 shutdown)

Case studies: Microbeam analysis

Resolution $3 \times 4 \mu\text{m}^2$

Scan amplitude $2.6 \times 2.6 \text{ mm}^2$

Operation modes: scan, raster, point



Detector configuration


Case studies: Microbeam analysis




ipfn
INSTITUTO DE PLASMAS
E FUSÃO NUCLEAR





W (M_{α})




W (L_{α})



Fe (K_{α})



$D(^3\text{He}, p_0)^4\text{He}$



$^9\text{Be}(^3\text{He}, p_n)^{11}\text{B}$

LP 26-15BN

W sphere

PIXE ~30min

NRA ~3h

