



Summary of high-resolution Littrow data analysis

Gennady Sergienko and Rongxing Yi
FZJ




This work has been carried out within the framework of the EUROfusion Consortium, partially funded by the European Union via the Euratom Research and Training Programme (Grant Agreement No 101052200 — EUROfusion). The Swiss contribution to this work has been funded by the Swiss State Secretariat for Education, Research and Innovation (SERI). Views and opinions expressed are however those of the author(s) only and do not necessarily reflect those of the European Union, the European Commission or SERI. Neither the European Union nor the European Commission nor SERI can be held responsible for them.

Littrow spectrometer for Balmer- α lines ratio measurement




Spectrometer set-up




Spectrometer parameters

- Focal length 750 mm
- Aperture ratio f/5.8
- Grating 1200 l/mm, 165x135 mm², blaze angle 57.22°
- Dispersion 0.685 nm/mm in 2nd order
- Slit 100 μ m (width) x 10 mm (height)
- Resolving power of 17900 at $\lambda = 656$ nm
- Line fiber bundle for slit side 46 fibers ($\varnothing 100$ μ m + 10 μ m (cladding))
- Etendue $2.4 \cdot 10^{-4}$ cm² srad
- Andor iStar CCD (DH334T-18F-A3) cooled intensified gated camera, 30% < QE < 40% at $\lambda = 600 - 700$ nm, 1024x1024 pix, 16bit, gating down to < 5 ns







Spectrometer calibration

- Wavelength calibration was performed with H_2/D_2 spectral lamp
- Absolute spectral sensitivity in 654-658 nm range was measured by calibrating integrating sphere placed in front of the LIBS tool cone
- Sensitivity curve produced by this calibration procedure is divided by calculated light collection efficiency of LIBS tool optics for 1 mm plasma size




Littrow spectral data






Highest D/H+D ratio

- Voigt fitting function (fixed Gaussian width, Lorentz width shared by H/D/T lines)
- Wavelength 6 pixels deviation
- No fixed initial value for H α position
- W I line interference make it difficult to find the accurate wavelength for H α

Fuel retention in module 2 divertor (poloidal from HFGC-Tile8)



Too weak D





Improvement of fitting formulae for Balmer- α line

Single Lorentz fit with shift:

$$f(\lambda) = y_0 + \frac{2A_0}{\pi} \frac{w_{FWHM}}{4(\lambda - \lambda_w)^2 + (w_{FWHM})^2}$$

$$\lambda_w = \lambda_0 + a (w_{FWHM})^{3/2} \quad \lambda_0 = 656.2714 \text{ nm} \quad a = 0.05877 \text{ nm}^{-1/2}$$

Double Lorentz fit with shift:

$$f(\lambda) = y_0 + \frac{2A_0}{\pi} \left(\frac{wC}{4(\lambda - \lambda_0 - aw^{3/2})^2 + w^2} + \frac{k w(1-C)}{4(\lambda - \lambda_0 - a(kw)^{3/2})^2 + (kw)^2} \right)$$

$$0 < y_0 \quad 0 < A_0 \quad 0 < w \quad 0 < k < 1 \quad 0 < C < 1$$

Tile 1 RH14W 2011-23 laser crater 1

